


## IV Labels


To configure the Pharmacy Plus to print IV Labels, you must have an IMPACT printer attached to LPT3 (a third parallel printer port). This printer will be used only for IV labels.

**Order labels:** IV Labels can be ordered via Rx Systems, Inc. (800) 922-9142 or Rep. Ed Tamara (954) 442-1628, ask for the Abacus IV labels (4 x 5).



To configure the impact printer on the Pharmacy Plus software, enter the Configuration Preference. Password ABACUS.

## Abacus – Pharmacy Plus: IV Labels



From the Configuration Preferences, press page down key (PgDn) until to reach the fifth page (above). Enter a printer for question C (Printer #3) and a form.

1. Printer: OKIDATA PRINTER
2. Form: VA

You are now ready to enter and print IV labels.

Abacus – Pharmacy Plus: IV Labels

Entering New/Refill prescriptions 3.98p

Patient#: ABDALL0000  
 Name: ABDALLAH, AIMAN  
 Address: [ Prescription entry NDC #49502-0697-60 ]  
 Ph: Drug: ALBUTEROL .83MG/ML S (DEY LABS.) Eq→PROVENTIL .83MG/ML S  
 Grp#: Pkg. Size: 3.00 AMP/Unit: 0.08 On-Hand: 0.00  
 Card: Doctor: ME0083251 BAAB, HOLLY ( )

HIPAA: Need Counseling  
 Generic: No

Plan	Rx No.	Written	Filled	P.Qty	D.Qty	Refills	tal	St
	NEW	08/10/2004	08/10/2004	1.000	1.000	0		

SIG code: INFUSE OVER 3 HOURS. Copies: 1 Origin: 0  
 Days supply: 1 Drug exp: 08/10/2005 Lot#:

IngrPrice	Disp.Fee	Discount	Total	CoPayment	Catg	Status
0.08	4.23	0.00	4.31	0.00		Not Billed

Line: 0.00

Use ' \ ' to concatenate- [AMP: 0.08 Acq: 0.06] Remaining Qty=0.000

F2 Est.\$ F3 DUR F4 ModDAT F5 IV F6 HCFA F7 WC F8 Page 2 F9 Notes PgDn Done

From the prescription entry screen, press the [F5] function key to get the IV screen.

Entering New/Refill prescriptions 3.98p

Patient#: ABDALL0000  
 Name: ABDALLAH, AIMAN  
 Address: [ IV LABEL ]  
 Ph: ENTER CODE OF IV LABEL TO RETRIEVE OR SAVE: IV  
 Grp#: Do MVI-CONC 5 ML VIAL 3.00 ML  
 Card: ASCORBIC ACID 25 GRAMS/50 ML 50.00 ML  
 Plan: Rx FOLINIC ACID 50 GRAMS/50 ML 5.00 MG  
 N VIT B12 (CYANOCOBALAMIN) 1000MCB 1.00 ML  
 SI NAEL 0.9% 500 ML 500.00 ML  
 In INFUSE OVER 3 HOURS.

HIPAA: Need Counseling  
 Generic: No

Plan	Rx No.	Written	Filled	P.Qty	D.Qty	Refills	tal	St

Line: 0.00

Cash Sale - [AMP: 1.98 Acq: 1.59] Remaining Qty=0.000

Esc Cancel F1 Help PgDn Done

From the IV screen you may enter a code to save the IV label information or you may skip the code field and begin to type what you want printed on the label. Once this step is done the rest is the same as a regular prescription.

Abacus – Pharmacy Plus: IV Labels

Below is a sample of an IV label (not to actual size).

XYZ Infusion Pharmacy 123 Sample Rd., Miami, FL 33111 (305) 555-1212	
Patient: ANDERSON, CHRISTINA Dr. ABAUNZA, RAMIRO J	
Rx 9902 ALBUTEROL .83MG/ML SOLUTION Inv #102584	
MVI-CONC 5 ML VIAL	3.00 ML
ASCORBIC ACID 25 GRAMS/50 ML	50.00 ML
FOLINIC ACID 50 MG VIAL	5.00 MG
VIT B12 (CYANDCOBALAMIN)	1.00 ML
NACL 0.9% 500 ML	500.00 ML
STORAGE: Refrigerate	
INFUSE OVER 3 HOURS.	
Prep Date 02/18/2000	Discard after: 02/01/2001
<b>Caution: Federal Law Prohibits Dispensing Without a Prescription</b>	